



Medical Microbiology

Third stage- College of Dentistry

Lab. 14

Neisseria spp.

Prepared by:

PhD. Dheyaa Naji Hamza

M.Sc. Zainab Majid Mohammed

Learning Objectives:

After completing this lab, you should be able to:

1-Describe the microscopic morphology and general Characteristics of the *Neisseria*.

4-Discuss the pathogenic classification and major patho-types of *Neisseria*.



2-Explain the laboratory techniques used for the detection and identification of *Neisseria*.

3-Identify the cultural characteristics of *Neisseria* and the appropriate media used for its isolation.



Morphology and General Characteristics of Neisseria:

- Shaped:** diplococci, typically described as **kidney-bean** or **coffee-bean** shaped with adjacent sides flattened.
- Gram stain:** Gram-negative bacteria.
- Respiration:** Aerobe and fastidious (requires enriched media) like **Chocolate Agar** for growth.
- Spore formation:** non-spore forming.
- Motility:** non- motile.
- Capsule:**
 - ☐ **Neisseria meningitidis:** capsulated, It has a polysaccharide capsule.
 - ☐ **Neisseria gonorrhoeae:** non-capsulated, though some strains may have a very thin capsule-like layer.



Habitat of *Neisseria* Species:

Neisseria species colonize the **mucous membranes of humans**.

-*Neisseria meningitidis*

Primarily inhabits the **nasopharynx (upper respiratory tract)** and may be present in healthy asymptomatic carriers.

-*Neisseria gonorrhoeae*

Mainly colonizes the **urogenital tract**.

-Commensal *Neisseria*

Some species are part of the **normal oral flora** and are generally non-pathogenic, including:

- *Neisseria sicca*
- *Neisseria mucosa*

These organisms typically do not cause disease under normal conditions.

Pathogenic *Neisseria* Species:

There are two major pathogenic species in humans:

1-*Neisseria gonorrhoeae* (Gonococcus):

Causes **gonorrhea**, Can also lead to **gonococcal pharyngitis**, which is important in dental practice due to oropharyngeal infection.

2-*Neisseria meningitidis* (Meningococcus)

Causes **meningitis** and **septicemia**.

Frequently found in the nasopharynx of asymptomatic carriers, facilitating transmission.



Virulence Factors:

1-Pili (Fimbriae):

Enable attachment to mucosal epithelial cells.

2-Capsule (mainly in *N. meningitidis*):


Polysaccharide capsule protects against phagocytosis.

3-Lipooligosaccharide (LOS):

Acts as endotoxin; triggers strong inflammatory response leading to tissue damage and septic shock.

4-IgA Protease:

Breaks down secretory IgA on mucosal surfaces, enhancing colonization.



Diagnosis Methods:


1-Culture Media:

They grow best at 37°C in an atmosphere enriched with 5-10% CO₂.

A-Blood Agar: Used mainly for *N. meningitidis*.

B-Chocolate Agar (CHOC): A non-selective enriched medium containing heated blood, which releases essential growth factors.

C-Modified Thayer-Martin (MTM) Agar: A selective medium based on **chocolate agar** with added **antibiotics** (Vancomycin, Colistin, Nystatin, and Trimethoprim) to inhibit the growth of Gram positive, Gram negative and fungi, allowing only *Neisseria* species to grow.



2-Biochemical Tests

1-Oxidase Test: positive.

2-Catalase Test: positive except *N. elongata*.

3-Carbohydrate Utilization Test (CTA Sugars):

Principle: Neisseria species are oxidative, not fermentative. They produce weak acids when they break down carbohydrates in the presence of oxygen.

Medium: CTA (Cystine Trypticase Agar) with 1% carbohydrate and Phenol Red as a pH indicator.

Results:

Positive: The medium turns Yellow (acid production) at the top of the tube.

Negative: The medium remains Red/Orange.

Species Differentiation:

- ***N. gonorrhoeae*:** Oxidizes **Glucose** only.
- ***N. meningitidis*:** Oxidizes **Glucose** and **Maltose**.
- ***N. lactamica*:** Oxidizes **Glucose**, **Maltose**, and **Lactose**.

N gonorrhoeae



N meningitidis



	Glucose	Maltose	Lactose
N. Gonorrhea	+	-	-
N. Meningitides	+	+	-



Any Question?



Thank you